

Human Chitinase 3-like 1/CHI3L1 protein (Recombinant) (His-Myc-Tag) (STJP018069)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	CHO Cells

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4)
Molecular Weight	43.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1116
Gene Symbol	CHI3L1
UniProt ID	CH3L1_HUMAN
Immunogen Region	22-383aa
Immunogen Sequence	YKLVCYYTSW SQYREGDGSC FPDALDRFLC THIIYSFANI SNDHIDTWEW NDVTLYGMLN TLKNRNPNLK TLLSVGGWNF GSQRFSKIAS NTQSRRTFIK SVPPFLRTHG FDGLDLAWLY PGRRDKQHFT TLIKEMKAEF IKEAQPGKKQ LLLSAALSAG KVTIDSSYDI AKISQHLDFI SIMTYDFHGA WRGTTGHHSP LFRGQEDASP DRFSNTDYAV GYMLRLGA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CHI3L1 protein 3µg under reducing conditions with Coomassie blue stain