

Human FXYD5/Dysadherin protein (Recombinant) (His-Tag) (STJP018075)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	16.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53827
Gene Symbol	FXYD5
UniProt ID	FXYD5_HUMAN
Immunogen Region	22-145aa
Immunogen Sequence	QTLKDTT SSSADSTIM DIQVPTRAPD AVYTELQPTS PTPTWPADET PQPQTQTQQL EGTGDLVTD PETHKSTKAA HPTDDTTLS ERPSSTDVQ TDPQTLKPSG FHEDDPFFYD EHTLRKR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human FXYD5 Dysadherin protein under reducing conditions