

Human delta Sarcoglycan/SGCD protein (Recombinant) (His-Tag) (STJP018087)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Proteins |
| Applications | SDS-PAGE/Denatured |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|---|
| Concentration | 1 mg/mL |
| Formulation | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol |
| Molecular Weight | 28.1 kDa |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | 6444 |
| Gene Symbol | SGCD |
| UniProt ID | SGCD_HUMAN |
| Immunogen Region | 57-289aa |
| Immunogen Sequence | KVMNFT IDGMGNLRIT EKGLKLEGDS EFLQPLYAKE IQSRPGNALY FKSARNVTVN ILNDQTKVLT QLITGPKAVE AYGKKFEVKT VSGKLLFSAD NNEVVVGAER LRVLGAEGTV FPKSIETPNV RADPFKELRL ESPTRSLVME APKGVEINAE AGNMEATCRT ELRLESKDGE IKLDAAKIRL PRLPHGSYTP TGTRQK |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human delta Sarcoglycan SGCD protein under reducing conditions