

## Human GRP/LGALSL protein (Recombinant) (His-Tag) (STJP018088)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	21.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">29094</a>
Gene Symbol	<a href="#">LGALSL</a>
UniProt ID	<a href="#">LEGL_HUMAN</a>
Immunogen Region	1-172aa
Immunogen Sequence	MAGSVAD SDAVVKLDDG HLNNLSLSSPV QADVFPRLI VPFCGHIKGG MRPGKKVLVM GIVDLNPESF AISLTCGDSE DPPADVAIEL KAVFTDRQLL RNCSISGERG EEQSAIPYFP FIPDQPFRVE ILCEHPRFRV FVDGHQLDFD YHRIQTL SAI DTIKINGDLQ ITKLG

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human GRP LGALSL protein under reducing conditions