

Human IL-34 protein (Recombinant) (His-Tag) (STJP018097)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	29.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146433
Gene Symbol	IL34
UniProt ID	IL34_HUMAN
Immunogen Region	21-242aa
Immunogen Sequence	NE PLEMWPLTQN EECTVTGFLR DKLQYRSRLQ YMKHYFPINY KISVPYEGVF RIANVTRLQR AQVSERELRY LWVLVLSAT ESVQDVLLEG HPSWKYLQEV QTLLLNQVQG LTDVEVSPKV ESVLSLLNAP GPNLKLVRPK ALLDNCFRVM ELLYCSCCKQ SSVLNWQDCE VPSPQSCSPE PSLQYA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human IL-34 protein SDS-PAGE gel under reducing conditions stained with Coomassie blue