

Human ID2 protein (Recombinant) (His-Tag) (STJP018098)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.4M urea
Molecular Weight	17 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3398
Gene Symbol	ID2
UniProt ID	ID2_HUMAN
Immunogen Region	1-134aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MKAFSPVRSV RKNLSLDHSL GISRSKTPVD DPMSLLYNMN DCYSKLKELV PSIPQNKKVS KMEILQHVID YIDLQIALD SHPTIVSLHH QRPQGNQASR TPLTTLNTDI SILSLQASEF PSELMSNDSK ALCG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3ug)

SDS-PAGE gel of recombinant human ID2 protein under reducing conditions stained with Coomassie blue