

Human CIAO1 protein (Recombinant) (His-Tag) (STJP018113)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	40 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9391
Gene Symbol	CIAO1
UniProt ID	CIAO1_HUMAN
Immunogen Region	1-339aa
Immunogen Sequence	MKDSLVL LGRVPAHPDS RCWFLAWNPA GTLLASCGGD RRIRIWGTEG DSWICKSVLS EGHQRTVRKV AWSPCGNYLA SASFDATTCT WKKNQDDFEC VTTLEGHENE VKSVAWAPSG NLLATCSRDK SVVWWEVDEE DEYECVSVLN SHTQDVKHVV WHPSQELLAS ASYDDTVKLY REEEDDWVCC ATLEGHESTV WSLAFD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CIAO1 protein under reducing conditions stained with Coomassie blue