

Human ARPC2 protein (Recombinant) (His-Tag) (STJP018117)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% Glycerol, 1mM DTT
Molecular Weight	36.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10109
Gene Symbol	ARPC2
UniProt ID	ARPC2_HUMAN
Immunogen Region	1-300aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSILLEVN NRIIEETLAL KFENAAAGNK PEAVEVTFAD FDGVLYHISN PNGDKTKVMV SISLKFKYKEL QAHGADELLK RVYGSFLVNP ESGYNVSLLY DLENLPASKD SIVHQAGMLK RNCFASVFEK YFQFQEEGKE GENRAVIHYR DDETMYESK KDRVTTVVFST VFKDDDDVVI GKVFMQEFKE GRRASHTAPQ VLFSHREP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ARPC2 protein 3 µg under reducing conditions
Coomassie stain