

Human CAMP/LL37/FALL39 protein (Recombinant) (His-Tag) (STJP018128)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% Glycerol
Molecular Weight	18.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	820
Gene Symbol	CAMP
UniProt ID	CAMP_HUMAN
Immunogen Region	34-173aa
Immunogen Sequence	QVLSYKE AVLRAIDGIN QRSSDANLYR LLDLDPRTM DGDPDTPKPV SFTVKETVCP RTTQQSPEDC DFKKDGLVKR CMGTVTLNQA RGSFDISCDK DNKR FALLGD FFRKSKEKIG KEFKRIVQRI KDFLRNLVPR TES

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CAMP/LL37 FALL39 protein under reducing conditions