

Human FAM49B protein (Recombinant) (His-Tag) (STJP018140)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.1M NaCl, 1mM DTT
Molecular Weight	39.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51571
Gene Symbol	CYRIB
UniProt ID	CYRIB_HUMAN
Immunogen Region	1-324aa
Immunogen Sequence	MGNLLKV LTCTDLEQGP NFFLDFENAIQ PTESEKEIYN QVNVVLKDAE GILEDLQSYR GAGHEIREAI QHPADEKLQE KAWGAVVPLV GKLLKFYEFS QRLEAALRGL LGALTSTPYS PTQHLEREQA LAKQFAEILH FTLRFDELKM TNPAIQNDFS YYRRTL SRMR INNVP AEGEN EVNNELANRM SLFYAEATPM LKTLSD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human FAM49B protein under reducing conditions with Coomassie blue stain