

Human CPOX protein (Recombinant) (His-Tag) (STJP018149)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	41.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1371
Gene Symbol	CPOX
UniProt ID	HEM6_HUMAN
Immunogen Region	111-454aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSTSLGRPE EEDELAHRC SSFMAPPVTD LGELRRRPGD MKTKMELLIL ETQAQVCQAL AQVDGGANFS VDRWERKEGG GGISCVLQDG CVFEKAGVSI SVVHGNLSEE AAQMRSRGK VLKTKDGKLP FCAMGVSSVI HPKNPHPTI HFNYRYFEVE EADGNKQWWF GGGCDLTPTY LNQEDAVHFFH RTLKEACDQH GPDLYPKF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CPOX protein 3µg under reducing conditions stained with Coomassie blue