

Human ID1 protein (Recombinant) (His-Tag) (STJP018153)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	18.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3397
Gene Symbol	ID1
UniProt ID	ID1_HUMAN
Immunogen Region	1-155aa
Immunogen Sequence	MKVASGS TATAAAGPSC ALKAGKTASG AGEVVRCLSE QSVASRCAG GAGARLPALL DEQQVNVLLY DMNGCYSLK ELVPTLPQNR KVSKEILQH VIDYIRDLQL ELNSESEVGT PGGRGLPVRA PLSTLNGEIS ALTAEAACVP ADDRILCR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ID1 protein
 3µg under reducing conditions Coomassie blue
 stain