

## Human RAB5C protein (Recombinant) (His-Tag) (STJP018154)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 1 mg/mL                                                                                                                                             |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% Glycerol                                                                         |
| Molecular Weight    | 25.9 kDa                                                                                                                                            |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                               |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">5878</a>                                                                                                                                                                                                          |
| Gene Symbol        | <a href="#">RAB5C</a>                                                                                                                                                                                                         |
| UniProt ID         | <a href="#">RAB5C_HUMAN</a>                                                                                                                                                                                                   |
| Immunogen Region   | 1-216aa                                                                                                                                                                                                                       |
| Immunogen Sequence | MAGRGA ARPNGPAAGN KICQFKLVLL GESAVGKSSL VLRVKGQFH EYQESTIGAA FLTQTVCLDD TTVKFEIWDI AGQERYHSLA PMYYRGAQAA IVVYDITNTD TFARAKNWWK ELQRQASPNI VIALAGNKAD LASKRAVEFQ EAQAYADDNS LLFMETSAKT AMNVNEIFMA IAKKLPKNEP QNATGAPGRN RGVDLQ |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RAB5C protein 3µg under reducing conditions Coomassie blue stain