

Human IGFLR1 protein (Recombinant) (His-Tag) (STJP018160)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.2M NaCl, 50% Glycerol, 2mM DTT
Molecular Weight	17.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79713
Gene Symbol	IGFLR1
UniProt ID	IGER1_HUMAN
Immunogen Region	23-163aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSSQYCGRL EYWNPDNKCC SSCLQRFGPP PCPDYEFREN CGLNDHGDV TPPFRKCSSG QCNPDAELC SPCGGGAVTP TPAAGGGRTT WRCRERPVPA KGHCPLTPGN PGAPSSQERS SPASSIAWRT PEPVPQQAWP NFLP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human IGFLR1 protein 3 µg reducing conditions Coomassie blue stain