

E.coli mug protein (Recombinant) (His-Tag) (STJP018164)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 20% Glycerol
Molecular Weight	21.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-168aa
Immunogen Sequence	MVEDILA PGLRVVFCGI NPGLSAGTG FPFAPANRF WKVIYQAGFT DRQLKPQEAQ HLLDYRCGVT KLVD RPTVQA NEVSKQELHA GGRKLIKIE DYQPQALAIL GKQAYEQGFS QRGAWGKQT LTIGSTQIWW LPNPSGLSRV SLEKLVEAYR ELDQALVVRG R

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E. coli mug protein under reducing conditions stained with Coomassie blue