

## Human CYB5R1 protein (Recombinant) (His-Tag) (STJP018173)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	33.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">51706</a>
Gene Symbol	<a href="#">CYB5R1</a>
UniProt ID	<a href="#">NB5R1_HUMAN</a>
Immunogen Region	29-305aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSLVRRSRR QVTLDPNE KYLLRLDKT TVSHNTKRFR FALPTAHTL GLPVGKHIYL STRIDGSLVI RPYTPVTSDE DQG YVDLVK VYKGVHPKF PEGGKMSQYL DSLKVGDVVE FRGPSGLLTY TGKGFNIQP NKKSPPEPRV AKKLGMIAGG TGITPMLQLI RAILKVPEDP TQCFLLFANQ TEKDILRED LEELQARY

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human CYB5R1 protein under reducing conditions stained with Coomassie blue