

Human ROBLD3/LAMTOR2 protein (Recombinant) (His-Tag) (STJP018174)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	15.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	28956
Gene Symbol	LAMTOR2
UniProt ID	LTIOR2_HUMAN
Immunogen Region	1-125aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMLRPKAL TQVLSQANTG GVQSTLLLNN EGSLLAYSGY GDTDARVTAA IASNIWAAYD RINGNQAFNE NLKFILMDCM EGRVAITRVA NLLLCMYAKE TVGFGMLKAK AQUALVQYLEE PLTQVAAS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



15% SDS-PAGE (3ug)

SDS-PAGE gel of recombinant human ROBLD3 LAMTOR2 protein under reducing conditions stained with Coomassie blue