

Human IL-17B protein (Recombinant) (His-Tag) (STJP018184)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	20 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27190
Gene Symbol	IL17B
UniProt ID	IL17B_HUMAN
Immunogen Region	21-180aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHMQPRSP KSKRKGQGRP GPLAPGPHQV PLDLVSRMKP YARMEEYERN IEEMVAQLRN SSELAQRKCE VNLQLWMSNK RSLSPWGYSI NHDPSPRPVD LPEARCLCLG CVNPFTMQED RSMVSPVFS QVPVRRRLCP PPPRTGPCRQ RAVMETIavg CTCIF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human IL-17B protein under reducing conditions stained with Coomassie blue