

Human RhoQ protein (Recombinant) (His-Tag) (STJP018188)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT
Molecular Weight	24.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23433
Gene Symbol	RHOQ
UniProt ID	RHOQ_HUMAN
Immunogen Region	1-202aa
Immunogen Sequence	MAHGPGA LMLKCVVVG D GAVGKTCLLM SYANDAFPEE YVPTVFDHYA VSVTVGGKQY LLGLYDTAGQ EDYDRLRPLS YPMTDVFLIC FSVVNPASFQ NVKEEWVPEL KEYAPNVVPL LIGTQIDLRD DPKTLARLND MKEKPICVEQ GQKLAK EIGA CCYVECSALT QKGLKTVFDE AIIAILTPKK HTVKKRIGSR CINCC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RhoQ protein 3µg under reducing conditions Coomassie stain