

## Human BMP-6 protein (Recombinant) (His-Tag) (STJP018197)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 10mM Sodium Citrate buffer (pH 3.5) containing 10% Glycerol
Molecular Weight	18 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">654</a>
Gene Symbol	<a href="#">BMP6</a>
UniProt ID	<a href="#">BMP6_HUMAN</a>
Immunogen Region	375-513aa
Immunogen Sequence	MSASSR RRQQRNRST QSQDVARVSS ASDYNSSELK TACRKHELYV SFQDLGWQDW IIAPKGYAAN YCDGECSPFL NAHMNATNHA IVQTLVHLMN PEYVPKPCCA PTKLNAISVL YFDDNSNVIL KKYRNMVVRA CGCH

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human BMP-6 protein under reducing conditions  
 Coomassie blue stain