

Human USE1 protein (Recombinant) (His-Tag) (STJP018206)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.1M NaCl, 1mM DTT
Molecular Weight	28.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55850
Gene Symbol	USE1
UniProt ID	USE1_HUMAN
Immunogen Region	1-231aa
Immunogen Sequence	MAASRL LNLVRLLSRC EAMAAEK RDP DEWRLEKYVG ALEDMLQALK VHASKPASEV INEYSWKVDF LKGM LQAEKL TSSSEKALAN QFLAPGRVPT TARERVPATK TVHLQSRARY TSEMRSELLG TDSAEP EMDV RKRTGVAGSQ PVSEKQSAE LDLVLQRHQN LQEKLA E EML GLARSLKTNT LAAQSVIKKD NQTL SH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human USE1 protein under reducing conditions with Coomassie blue stain