

Human RPL23A protein (Recombinant) (His-Tag) (STJP018209)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT, 250mM imidazole
Molecular Weight	20.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6147
Gene Symbol	RPL23A
UniProt ID	RL23A_HUMAN
Immunogen Region	1-156aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSMAPKAKK EAPAPPKAEA KAKALKAKKA VLKGVHSHKK KKIRTSPFTR RPKTLRLRRQ PKYPRKSAPR RNKLDHYAII KFPLTTESAM KKIEDNNTLV FIVDVKANKH QIKQAVKKLY DIDVAKVNTL IRPDGEKKAY VRLAPDYDAL DVANKIGII

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPL23A protein 3ug under reducing conditions