

Human ABRACL protein (Recombinant) (His-Tag) (STJP018210)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol
Molecular Weight	11.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58527
Gene Symbol	ABRACL
UniProt ID	ABRAL_HUMAN
Immunogen Region	1-81aa
Immunogen Sequence	MNVDHEV NLLVEEIHRL GSKNADGKLS VKFGVLFRRDD KCANLFEALV GTLKAARRK IVTYPGELL QGVHDDVDII LLQD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ABRACL protein 3µg under reducing conditions Coomassie stain