

Human RABL5/IFT22 protein (Recombinant) (His-Tag) (STJP018216)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	23.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64792
Gene Symbol	IFT22
UniProt ID	IFT22_HUMAN
Immunogen Region	1-185aa
Immunogen Sequence	MLKAKIL FVGPCESGKT VLANFLTESS DITEYSPTQG VRILEFENPH VTSNNKGTGC EFELWDCGGD AKFESCWPAL MKDAHGVVIV FNADIPSHRK EMEMWYSCFV QPQLQDTQC MLIHHKPGS GDDKGSLSLS PPLNKLKLVH SNLEDDPEEI RMEFIKYLKS IINSMESRD REEMSIMT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RABL5 IFT22 protein under reducing conditions
Coomassie blue stain