

Human ACOT8 protein (Recombinant) (His-Tag) (STJP018222)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 2mM DTT
Molecular Weight	38.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10005
Gene Symbol	ACOT8
UniProt ID	ACOT8_HUMAN
Immunogen Region	1-319aa
Immunogen Sequence	MSSPQAP EDGQGCGRG DPPGDLRSVL VTTVLNLEPL DEDLFRGRHY WVPKRLFGG QIVGQALVAA AKSVSESVHV HSLHCYFVRA GDPKLPVLYQ VERTRTGSSF SVRSVKAVQH GKPIFCQAS FQQAQPSMQ HQFSMPTVPP PEELLDCELT IDQYLDPNL QKRYPLALNR IAAQEVPIE KPVNPSPLSQ LQRMPEP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human ACOT8 protein under reducing conditions Coomassie blue stain