

Human FAM50A protein (Recombinant) (His-Tag) (STJP018226)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	25.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9130
Gene Symbol	FAM50A
UniProt ID	FA50A_HUMAN
Immunogen Region	150-339aa
Immunogen Sequence	TTKKRKL GKNPDVDTSF LPDRDREEEE NRLREELRQE WEAKQEIKIS EEIEITFSYW DGSGHRRTVK MRKGNTMQQF LQKALEILRK DFSELSAGV EQLMYIKEDL IIPHHHSFYD FIVTKARGKS GPLFNFDVHD DVRLLSDATV EKDESHAGKV VLRSWYEKKNK HIFPASRWEP YDPEKKWDKY TIR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FAM50A protein 3µg under reducing conditions stained with Coomassie blue