

Human CCDC90B protein (Recombinant) (His-Tag) (STJP018227)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	24.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60492
Gene Symbol	CCDC90B
UniProt ID	CC90B_HUMAN
Immunogen Region	43-230aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MSGYDRRPV DITPLEQRKL TFDTHALVQD LETHGFDKTQ AETIVSALTA LSNVSLDTIY KEMVTQAQQE ITVQQLMAHL DAIRKDMVIL EKSEFANLRA ENEKMKIELD QVKQQLMHET SRIRADNKLD INLERSRVT D MFTDQEKQLM ETTTEFTKKD TQTKSISET SNKIDAEIAS LKTLMESNKL E

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel showing recombinant human CCDC90B protein 3µg under reducing conditions stained with Coomassie blue