

Human FIBP protein (Recombinant) (His-Tag) (STJP018252)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	44.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9158
Gene Symbol	FIBP
UniProt ID	FIBP_HUMAN
Immunogen Region	1-364aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMTSELDI FVGNTTLIDE DVYRLWLDGY SVTDAVALRV RSGILEQTGA TAAVLQSDTM DHYRTFHMLE RLLHAPPKLL HQLIFQIPPS ROALLIERYY AFDEAFVREV LGKKLSKGTK KDLDISTKT GITLKSCRRQ FDNFKRVFKV VEEMRGLVD NIQQHFLSD RLARDYAAIV FFANNRFETG KKKLQYLSFG DFAFCAEL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FIBP protein under reducing conditions Coomassie blue stain