

Human C1D protein (Recombinant) (His-Tag) (STJP018254)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT
Molecular Weight	18.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10438
Gene Symbol	C1D
UniProt ID	C1D_HUMAN
Immunogen Region	1-141aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSMAGEEIN EDYPVEIHEY LSAFENSIGA VDEMLKTMMMS VSRNELLQKL DPLEQAKVDL VSAYTLNSMF WYVLATQGVN PKEHPVKQEL ERIRVYMN RV KEITDKKKAG KLDRGAASRF VKNALWEPKS KNASKVANKG KSKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human C1D protein 3ug under reducing conditions with Coomassie blue stain