

Human CPPED1 protein (Recombinant) (His-Tag) (STJP018255)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol
Molecular Weight	37.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55313
Gene Symbol	CPPED1
UniProt ID	CPPED_HUMAN
Immunogen Region	1-314aa
Immunogen Sequence	MSAAEAG GVFHRARGRT LAAPPAEKES EWKGPFYFIL GADPQFGLIK AWSTGDCDNG GDEWEQEIRL TEQAVQAINK LNP KPFFVL CGDLIHAMPK KPWRTEQTED LKRVLRVDR AIPLVLVSGN HDIGNTPAE TVEEFCRTWG DDYFSFWVGG VLFLVLSQF YENPSKCPSL KQAQDQWLDE QLSIARQRHC QHAIVF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CPPED1 protein 3µg under reducing conditions Coomassie blue stain