

Human HBG1 protein (Recombinant) (His-Tag) (STJP018260)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol
Molecular Weight	18 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3047
Gene Symbol	HBG1
UniProt ID	HBG1_HUMAN
Immunogen Region	1-147aa
Immunogen Sequence	MGHFTEE DKATITSLWG KVNVEDAGGE TLGRLLVVYP WTQRRFFDSFG NLSSASAIMG NPKVKAHGKK VLTSLGDATEK HLDDLKGTFA QLSELHCDKL HVDPENFKLL GNVLVTVLAI HFGKEFTPEV QASWQKMTA VASALSSRYH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HBG1 protein 3µg under reducing conditions stained with Coomassie blue