

Human COX5B protein (Recombinant) (His-Tag) (STJP018263)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	13 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1329
Gene Symbol	COX5B
UniProt ID	COX5B_HUMAN
Immunogen Region	32-129aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSASGGGVP TDEEQATGLE REIMLAACKG LDPYNVLAPK GASGTREDPN LVPSISNKRI VGCICEEDNT SVVWFHLHKG EAQRCPRCGA HYKLVPQQLA H

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human COX5B protein under reducing conditions
 Coomassie blue stain