

Human RAB8B protein (Recombinant) (His-Tag) (STJP018282)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	25.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51762
Gene Symbol	RAB8B
UniProt ID	RAB8B_HUMAN
Immunogen Region	1-204aa
Immunogen Sequence	MAKTYDY LFKLLIGDS GVGKTCLLFR FSEDAFNITF ISTIGIDFKI RTIELDGKKI KLQIWDTAGQ ERFRTITTAY YRGAMGIMLV YDITNEKSFN NIKNWIRNIE EHASSDVERM ILGNKCDMND KRQVSKERGE KLAIDYGIKF LETSAKSSAN VEEAFFTLAR DIMTKLNRKM NDSNSAGAGG PVKITENRSK KTSFFR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant human RAB8B protein SDS-PAGE
gel 3 µg Coomassie blue stain