

Human BEX1 protein (Recombinant) (His-Tag) (STJP018284)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% Glycerol
Molecular Weight	17.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55859
Gene Symbol	BEX1
UniProt ID	BEX1_HUMAN
Immunogen Region	1-125aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSMSKEKR AVNSLSMENA NQENEEKEQV ANKGEPLALP LDAGEYCVPR GNRFRFRVRQ PILQYRWDMM HRLGEPQARM REENMERIGE EVRQLMEKLR EKQLSHSLRA VSTDPHHDH HDEFCLMP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human BEX1 protein under reducing conditions stained with Coomassie blue