

## Human Prohibitin 2 protein (Recombinant) (His-Tag) (STJP018290)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% Glycerol
Molecular Weight	35.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">11331</a>
Gene Symbol	<a href="#">PHB2</a>
UniProt ID	<a href="#">PHB2_HUMAN</a>
Immunogen Region	1-299aa
Immunogen Sequence	MAQNLKD LAGRLPAGPR GMGTALKLLL GAGAVAYGVR ESVFTVEGGH RAIFFNRIGG VQQDTILAEG LHFRIPIWFQY PIIYDIRARP RKISSPTGSK DLQMVNISR VLSRPNAQEL PSMYQRLGLD YEERVLPSIV NEVLKSVVAK FNASQLITQR AQVSLIRRE LTERAKDFSL ILDDVAITEL SFSREYTA AV EAKQVA

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Prohibitin 2 protein under reducing conditions stained with Coomassie blue