

Human HLA-DOA protein (Recombinant) (His-Tag) (STJP018292)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% Glycerol
Molecular Weight	24.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3111
Gene Symbol	HLA-DOA
UniProt ID	DOA_HUMAN
Immunogen Region	26-217aa
Immunogen Sequence	TKADHMG SYGPAFYQSY GASGQFTHF DEEQLFSVDL KKSEAVWRP EFGDFARFDP QGGLAGIAAI KAHLDILVER SNRSRAINVP PRVTVLPKSR VELGQPNILI CIVDNIFFPV INITWLRNGQ TVTEGVAQTS FYSQPDHLFR KFHLYLPFVPS AEDVYDCQVE HWGLDAPLLR HWELQVPIPP PDAME

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HLA-DOA protein 3ug reducing conditions Coomassie blue stain