

## Human HNRNPK protein (Recombinant) (His-Tag) (STJP018296)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 1 mg/mL                                                                                                                                             |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT                                                                |
| Molecular Weight    | 33 kDa                                                                                                                                              |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                      |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">3190</a>                                                                                                                                                                                                                                                 |
| Gene Symbol        | <a href="#">HNRNPK</a>                                                                                                                                                                                                                                               |
| UniProt ID         | <a href="#">HNRPK_HUMAN</a>                                                                                                                                                                                                                                          |
| Immunogen Region   | 1-276aa                                                                                                                                                                                                                                                              |
| Immunogen Sequence | MGSSHHHHHH SSGLVPRGSH MGSMETEQPE ETFPNTETNG EFGKRPAEDM EEEQAFKRSP<br>NTDEMVELRI LLQSKNAGAV IGKGGKNIKA LR TDYNASVS VPDSSGPERI LSISADIETI GEILKKIIP<br>LEEGLQLPSP TATSQLPLES DAVECLNYQH YKGSDFDCEL RLLIHQSLAG GIIGVKGAKI KELRENTQTT<br>IKLFQECCPH STD RVVLIGG KPDRVVEC |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HNRNPK protein under reducing conditions Coomassie blue stain