

Human MEOX2 protein (Recombinant) (His-Tag) (STJP018300)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 2mM DTT
Molecular Weight	15.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4223
Gene Symbol	MEOX2
UniProt ID	MEOX2_HUMAN
Immunogen Region	188-304aa
Immunogen Sequence	RKERTAF TKEQIRELEA EFAHHNYLTR LRRYEIAVNL DLTERQVKVW FQNRMRKWK R VKGGQQGAAA REKELVNVKK GTLLPSEL SG IGAATLQQTG DSIANEDSHD SDHSSEHAHL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human MEOX2 protein under reducing conditions Coomassie blue stain