

Human RPL22 protein (Recombinant) (His-Tag) (STJP018308)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	17 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6146
Gene Symbol	RPL22
UniProt ID	RL22_HUMAN
Immunogen Region	1-128aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAPVKKL VVKGKKKKKQ VLKFTLDCTH PVEDGIMDAA NFEQFLQERI KVNGKAGNLG GGVVTIERSK SKITVTSEVP FSKRYLKYL KKYLKNNLR DWLRVVANSK ESYELRYFQI NQDEEEEEDE D

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPL22 protein under reducing conditions stained with Coomassie blue