

Human RPS16 protein (Recombinant) (His-Tag) (STJP018316)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	18.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6217
Gene Symbol	RPS16
UniProt ID	RS16_HUMAN
Immunogen Region	1-146aa
Immunogen Sequence	MPSKGPL QSVQVFGRKK TATAVAHCKR GNGLIKVNGR PLEMIEPRTL QYKLLPEVLL LGKERFAGVD IRVRVKGGGH VAQIYAIRQS ISKALVAYYQ KYVDEASKKE IKDILIQYDR TLLVADPRRC ESKKFGGPGA RARYQKSYR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPS16 protein 3µg under reducing conditions stained with Coomassie blue