

Human IL-17F protein (Recombinant) (His-Tag) (STJP018320)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing, 10% Glycerol, 0.4M urea
Molecular Weight	17.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	112744
Gene Symbol	IL17F
UniProt ID	IL17F_HUMAN
Immunogen Region	31-163aa
Immunogen Sequence	RKIPK VGHTFFQKPE SCPPVPGGSM KLDIGIINEN QRVMSRNIE SRSTSPWNYT VTWDPNRYPS EVVQAQCRNL GCINAQGGKED ISMNSVPIQQ ETLVVRKHQ GCSVSFQLEK VLVTVGCTCV TPVIHHVQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human IL-17F protein under reducing conditions