

Human RPL11 protein (Recombinant) (His-Tag) (STJP018325)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 1mM DTT, 2mM EDTA, 250mM Imidazole
Molecular Weight	22.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6135
Gene Symbol	RPL11
UniProt ID	RL11_HUMAN
Immunogen Region	1-178aa
Immunogen Sequence	MAQDQGE KENPMRELRI RKLCLNICVG ESGDRLTRAA KVLEQLTGQT PVFSKARYTV RSFGIRRNEK IAVHCTVRGA KAEIILEKGL KVREYELRKN NFSDTGNFGF GIQEHIDLGI KYDPSIGIYG LDFYVVLGRP GFSIADKKRR TGCIGAKHRI SKEEAMRWFG QKYDGIILPG K

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing 3ug recombinant human RPL11 protein under reducing conditions. Coomassie blue stained.