

Human EMC2 protein (Recombinant) (His-Tag) (STJP018332)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 40% Glycerol, 2mM DTT
Molecular Weight	37.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9694
Gene Symbol	EMC2
UniProt ID	EMC2_HUMAN
Immunogen Region	1-297aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMKVSSEL YDVTWEEMRD KMRKWRENS RNSEQIVEVG EELINEYASK LGDDIWIIE QVMIAALDYG RDDLALFCLQ ELRRQFPQSH RVKRLTGMRP EAMERYDDAI QLYDRILQED PTNTAARKRK IAIRKAQGNK VEAIRELNEY LEQFVGQEA WHELAELYIN EHDYAKAFC LEELMMTNPH NHLYCQYAE VKYTQGG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human EMC2 protein under reducing conditions stained with Coomassie blue