

Human LIMD2 protein (Recombinant) (His-Tag) (STJP018339)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol
Molecular Weight	16.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80774
Gene Symbol	LIMD2
UniProt ID	LIMD2_HUMAN
Immunogen Region	1-127aa
Immunogen Sequence	MFQAAGA AQATPSHDAK GGGSSTVQRS KSFSRLRAQVK ETCAACQKTV YPMERLVADK LIFHNSCFCC KHCHTKLSLG SYAALHGEFY CKPHFQQLFK SKGNYDEGFG RKQHKELWAH KEVDPGTKTA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human LIMD2 protein 3ug under reducing conditions
 Coomassie blue stain