

Human NCK2 protein (Recombinant) (His-Tag) (STJP018340)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	45.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8440
Gene Symbol	NCK2
UniProt ID	NCK2_HUMAN
Immunogen Region	1-380aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMTTEEVIV IAKWDYTAQQ DQELDIKKNE RLWLLDDSKT WWRVRNAANR TGYVPSNYVE RKNSLKKGSL VKNLKDTLGL GKTRRKTSAR DASPTPSTDA EYPANGSGAD RIYDLNIPAF VKFAYVAERE DELSLVKGSR VTMKEKCS DG WWRGSYNGQI GWFPSNYVLE EVDEAAAE SP SFLSLRKGAS LSNGQGSRVL HVVQTLYP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human NCK2 protein 3µg under reducing conditions
Coomassie blue stain