

Human HOXC11 protein (Recombinant) (His-Tag) (STJP018376)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% Glycerol
Molecular Weight	36.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3227
Gene Symbol	HOXC11
UniProt ID	HXC11_HUMAN
Immunogen Region	1-304aa
Immunogen Sequence	MFNSVNL GNFCSPSRKE RGADFGERS CASNLYLPSC TYMPEFSTV SSFLPQAPSR QISYPYSAQV PPVREVS YGL EPSGKWHHRN SYSSCYAAAD ELMHRECLPP STVTEILMKN EGSYGGHHHP SAPHATPAGF YSSVNKNSVL PQAfDRFFDN AYCgGGDPPA EPPCSGKGEA KGEPEAPPAS GLASRAEAGA EAEAE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant human HOXC11 protein SDS-PAGE gel 3µg reducing conditions Coomassie blue stain