

Human TP53AIP1 protein (Recombinant) (His-Tag) (STJP018382)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	15kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	63970
Gene Symbol	TP53AIP1
UniProt ID	TPIP1_HUMAN
Immunogen Region	1-124aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSSSEASFR SAQASC SGAR RQGLGRGDQN LSVMP PNGRA QHTHPGWVSD PLVLGAQVHG GCRGIEALSV SSGSWSSATV WILTGLGLGL SRPFLPGATV LRDRPLGSAF ELSYDQKKAP LRLQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human TP53AIP1 protein under reducing conditions stained with Coomassie blue