

Human RNF7 protein (Recombinant) (His-Tag) (STJP018390)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	15.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9616
Gene Symbol	RNF7
UniProt ID	RBX2_HUMAN
Immunogen Region	1-113aa
Immunogen Sequence	MADVEDG EETCALASHS GSSGSKSGGD KMFSLKKWNA VAMWSWDVEC DTCAICRVQV MDACLRCQAE NKQEDCVVVW GECNHSFHNC CMLWVKQNN RCPLCQQDWV VQRIGK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human RNF7 protein under reducing conditions stained with Coomassie blue