

## Human NEURL2 protein (Recombinant) (His-Tag) (STJP018400)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Proteins           |
| Applications  | SDS-PAGE/Denatured |
| Host / Source | E.coli             |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.5 mg/mL                                                                                                                                          |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% Glycerol                                                                         |
| Molecular Weight    | 34.1 kDa                                                                                                                                           |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">140825</a>                                                                                                                                                                                                         |
| Gene Symbol        | <a href="#">NEURL2</a>                                                                                                                                                                                                         |
| UniProt ID         | <a href="#">NEUL2_HUMAN</a>                                                                                                                                                                                                    |
| Immunogen Region   | 1-285aa                                                                                                                                                                                                                        |
| Immunogen Sequence | MAAASEP VDSGALWGLE RPEPPPTRFH RVHGANIRVD PSGTRATRVE SFAHGVCFSR EPLAPGQVFL VEIEEKELGW CGHLRLGLTA LDPASLAPVP EFSLPDLVNL GHTWVFAITR HHNRVPREGR PEAEAAAPSR PPTLLVEPYL RIEQFRIPRD RLVGRSRPGL YSHLLDQLYE LNVLPPTARR SRLGVLCPR PDGTAD |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NEURL2 protein 3µg under reducing conditions Coomassie stain