

Human CUEDC2 protein (Recombinant) (His-Tag) (STJP018402)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	34.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79004
Gene Symbol	CUEDC2
UniProt ID	CUED2_HUMAN
Immunogen Region	1-287aa
Immunogen Sequence	MELERIV SAALLAFVQT HLEADLSGL DEVIFSYYLG VLEDLGPSGP SEENFDMEAF TEMMEAYVPG FAHIPRGITG DMMQKLSGQL SDARNKENLQ PQSSGVQGVV PISPEPLQRP EMLKEETRSS AAAAADTQDE ATGAEEELLP GVDVLLVFP TCSVEQAQWV LAKARGDLEE AVQMLVEGKE EGPAAWEGPN QDLPRR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human CUEDC2 protein under reducing conditions with Coomassie blue stain